

***Information Sheet. Product Code:4239/G Product: Lyndhurst Folding Shower Seat***



**User weight  
limit  
160kg  
(25 stones)**

Should you have a problem with this device, please quote the serial number attached to the frame

**DESCRIPTION:**

A permanently attached shower seat with padded back and arms, which can be lifted back against the wall, when not in use. The padded plastic seat provides extra comfort and can be removed for cleaning. Support legs are provided, with adjustment to suit the various heights of the users.

**INTENDED PURPOSE:**

As description. Will allow the user to sit in comfort and security whilst taking a shower. The 4239G has a cut out to allow easier bathing.

**INSTRUCTIONS FOR USE:**

To lower, pull on the top edge of the seat downwards then lower each arm. The legs will automatically fold outwards into position. After use, lift the seat and arms back into the vertical position and the device will hold in place.

The seat may be cleaned with domestic cleaning agents and detergents. Do not use abrasive materials as this may damage the protective frame coating.

***DO NOT lower the seat by pulling on the legs, as this will damage the stainless steel link rod and the leg itself. Warranty will be invalidated if this instruction is not adhered to.***

**FITTING INSTRUCTIONS:**

This device must be fitted to a load-bearing wall by a competent person and fastened with expansion bolts. It is essential that the legs are adjusted for height before the seat is fastened to the wall. After fitting, the legs must be vertical. The legs may be adjusted individually by unscrewing the threaded leg and locking it into position by tightening the nut.

**SAFETY WARNING:**

**Do not stand on this device.**

**Ensure the user is competent to use the product.**